

Five Fingered Disposable Nitrile Gloves EN 455 & EN 374 / CE / Non-Sterile / Powder Free / Category III

ILB-Product ID-No.: GL SYK Nitrile-np-non



Item	Details	Remark
CE Certificate	SANTRA	NB 2777
Classification	CE acc. to PPE Regulation EU 2016-425	Category III
Classification	CE acc. to Directive 93/42/EEC	Class I
Test Report by	SGS and Pony Testing	Swiss Testing Institution
Applied standard	EN 374, EN 455	
Factory ISO Certificates	ISO 9001, ISO 13485	
Product classification	Medical Use	
Sterile / Non Sterile	Non Sterile	

Color	Blue	Different size available
Packaging	1000 pcs/ Ctn 4.5 kg / Ctn	

Five Fingered Disposable Nitrile Gloves

This product is made of high quality nitrile material, not to tear easily, comfortable and smooth touch in hand, easy to wear, no holes and leakage. The CE certificate is issued by SANTRA, with notified body number of 2777, and the test report is done by SGS and Pony Testing

TECHNICAL SPECIFICATIONS

- CE 2777
- Gloves open and slide on easily, allowing for greater convenience and ease of use
- The shape ensures that you still have maximum range of movement with the gloves on
- Gloves extend to cover part of the wrist as well to provide further protection from splashes and messes
- Product size: XS/S/M/L/XL
- Material: high quality nitrile
- Packing unit: pack of 1000 pcs

Specifications:

1. EN 455-2-2015 Medical gloves for single use-part 2

Testing for physical properties: pass

2. EN 374-1: 2003 Minimum liquid proof length: Pass

3. EN 374-2:2014 protective gloves against chemicals and microorganisms: determination of resistance penetration: Pass

4. EN 374-3: 2013 Resistance to permeation by chemicals: Pass

Citric Acid: Class 6

Sodium hydroxide: Class 6

Isopropyl alcohol: class 0

Hexane: class 0

5. EN 374-4: 2013 Resistance by degradation by chemicals: Pass

6. EN 374-5:2016 Protection against bacteria and fungi, virus: pass