

Disposable Isolation Gown EN 13795-2 / Class I / CE / Non-Sterile / No Tape

ILB-Product ID-No.: PS SUNG-non-no-tape



Item	Details	Remark
China Ministry of Commerce White List Producer	Yes	CE
CE Certificate	Lotus NL BV as EU representative	EC-REP
Test Report by	Intertek	
Applied standard	EN 13795-2	
Factory ISO Certificates	ISO 9001 / 14001	
Product classification	Medical Use	
Sterile / Non Sterile	Non Sterile	
Color	Blue	Different size available
Packaging	100 pcs/ ctn 15 kg / ctn	

Disposable Isolation Gown, non-sterile, no tape

This Isolation Gown provide barrier thus protection for the personnel in non-medical situation with potential risk of contacting infectious patients, or particulates in the air, etc. It is also applicable to the general isolating of medical institutions in outpatient service area, inpatient ward, inspection room. Lotus NL BV as EU representative, and the test report is done by Intertek.

TECHNICAL SPECIFICATIONS

- CE, EC-REP
- Constructed from PE coated Non-woven fabric
- Offers economical, comfortable and reliable protection
- Sewn seams provide superior strength
- Product size: M/XL/3XL
- This non taped, non-sterile disposable surgical gown used for hospital surgery, treatment of patients and home care. It manufactured by high-quality PE coated non-woven fabrics. So, it has features of comfortable, soft, and non-irritant. This can protect against infection effective. It is the ideal choice for patient surgical and fluid isolation in hospitals and clinics.
- Material: PE coated Non-woven fabric, 41g
- Packing unit: pack of 100 pcs